

The Accutech table is built from sturdy 14 gauge tube construction and available in knocked down or fully welded construction. With it's height adjustability, numerous standard sizes, and custom sizing available, it's the most versatile table system on the market. The Accutech table gives your lab or educational space the ability to adapt to changing needs while easily pairing with our Infinity Bench and SteelTech cabinets.



PERFORMANCE AND AFFORDABILITY - WITHOUT COMPROMISE

- Chemical resistant epoxy powder-coat
- Built with 14-gauge precision laser cut tubing
- Standard widths of 48", 60", 72", 84" and 96"
- 29"-36" height adjustable
- Custom sizes upon request

OPTIONAL EQUIPMENT

- Modesty panel
- Modesty panel with wire organizer
- Choice of leveling glides or casters
- Electrical receptacles
- Data receptacles
- Knocked down or fully welded construction
- Steel, epoxy, phenolic and laminate top options
- Optional lower leg stringers for added rigidity

