

SELF-CLOSING FLAMMABLE STORAGE CABINETS

Newly redesigned for 2021, flammable storage cabinets are made in accordance with OSHA and NFPA standards for the storage of flammable and combustible liquids. For maximum protection of personnel and products, self-closing doors shut and latch automatically in the event of a fire.



Standard Features

- Three-point, self-latching doors
- Recessed paddle handle
- Self-closing door hinges (for obstruction free-access to contents)
- Adjustable steel shelf and leak-proof drip pan
- Heavy-duty, 18-gauge steel
- Double-walled construction
- Insulated cabinet frame and doors
- Upper and lower vents with threaded metal plugs
- Adjustable leveling glides
- Large “FLAMMABLE” warning label
- Durable powder coat finish with excellent chemical resistance



In the event of a fire, self-close doors shut and latch automatically when the fusible link melts at 165°F.



Part Numbers – 18” Deep

Part Number	Height	Width	Depth	Doors	Hinge Side	Weight (lbs)	
-------------	--------	-------	-------	-------	------------	--------------	--

29” High (Seated Height) x 18” Deep Self-Closing Flammable Storage Cabinets

FLAM2918L-18	29	18	18	Single	Left	100	
FLAM2918R-18	29	18	18	Single	Right	100	
FLAM2924L-18	29	24	18	Single	Left	110	
FLAM2924R-18	29	24	18	Single	Right	110	
FLAM2930-18	29	30	18	Double	Both	160	
FLAM2936-18	29	36	18	Double	Both	170	
FLAM2942-18	29	42	18	Double	Both	180	
FLAM2948-18	29	48	18	Double	Both	190	

32.5” High (ADA Height) x 18” Deep Self-Closing Flammable Storage Cabinets

FLAM3318L-18	32.5	18	18	Single	Left	110	
FLAM3318R-18	32.5	18	18	Single	Right	110	
FLAM3324L-18	32.5	24	18	Single	Left	120	
FLAM3324R-18	32.5	24	18	Single	Right	120	
FLAM3330-18	32.5	30	18	Double	Both	170	
FLAM3336-18	32.5	36	18	Double	Both	180	
FLAM3342-18	32.5	42	18	Double	Both	190	
FLAM3348-18	32.5	48	18	Double	Both	200	

35” High (Standing Height) x 18” Deep Self-Closing Flammable Storage Cabinets

FLAM3518L-18	35	18	18	Single	Left	110	
FLAM3518R-18	35	18	18	Single	Right	110	
FLAM3524L-18	35	24	18	Single	Left	120	
FLAM3524R-18	35	24	18	Single	Right	120	
FLAM3530-18	35	30	18	Double	Both	170	
FLAM3536-18	35	36	18	Double	Both	180	
FLAM3542-18	35	42	18	Double	Both	190	
FLAM3548-18	35	48	18	Double	Both	200	



Single Door, Hinge on Right



Double Door

Part Numbers – 22” Deep

Part Number	Height	Width	Depth	Doors	Hinge Side	Weight (lbs)	
-------------	--------	-------	-------	-------	------------	--------------	--

29” High (Seated Height) x 22” Deep Self-Closing Flammable Storage Cabinets

FLAM2918L-22	29	18	22	Single	Left	105	
FLAM2918R-22	29	18	22	Single	Right	105	
FLAM2924L-22	29	24	22	Single	Left	115	
FLAM2924R-22	29	24	22	Single	Right	115	
FLAM2930-22	29	30	22	Double	Both	165	
FLAM2936-22	29	36	22	Double	Both	175	
FLAM2942-22	29	42	22	Double	Both	185	
FLAM2948-22	29	48	22	Double	Both	195	

32.5” High (ADA Height) x 22” Deep Self-Closing Flammable Storage Cabinets

FLAM3318L-22	32.5	18	22	Single	Left	115	
FLAM3318R-22	32.5	18	22	Single	Right	115	
FLAM3324L-22	32.5	24	22	Single	Left	125	
FLAM3324R-22	32.5	24	22	Single	Right	125	
FLAM3330-22	32.5	30	22	Double	Both	175	
FLAM3336-22	32.5	36	22	Double	Both	185	
FLAM3342-22	32.5	42	22	Double	Both	195	
FLAM3348-22	32.5	48	22	Double	Both	205	

35” High (Standing Height) x 22” Deep Self-Closing Flammable Storage Cabinets

FLAM3518L-22	35	18	22	Single	Left	115	
FLAM3518R-22	35	18	22	Single	Right	115	
FLAM3524L-22	35	24	22	Single	Left	125	
FLAM3524R-22	35	24	22	Single	Right	125	
FLAM3530-22	35	30	22	Double	Both	175	
FLAM3536-22	35	36	22	Double	Both	185	
FLAM3542-22	35	42	22	Double	Both	195	
FLAM3548-22	35	48	22	Double	Both	205	



Single Door, Hinge on Right



Double Door