



BASIC AIRFLOW DIAGRAM—BV SERIES OPEN BYPASS

OPEN BYPASS PROVIDES PRIMARY AIRFLOW FROM ABOVE THE HOOD WITH SASH CLOSED

WITH SASH IN OPEN POSITION BYPASS AIRFLOW IS RESTRICTED AND FACE OF HOOD PROVIDES PRIMARY AIRFLOW

RAISED AIRFOIL PROVIDES AIRFLOW AT FACE OF HOOD WHEN SASH IS CLOSED

