

To meet the growing demand for flexibility in the laboratory furniture space, iQ Labs has developed our Infinity Bench system. With integrated plumbing and electrical service, quick disconnect options, and shared frame design, the height adjustable iQ Infinity Bench is the new standard for laboratory furniture from iQ Labs. Built around the need to quickly adjust and adapt to changes in the laboratory space the Infinity Bench is meant to provide ease of installation, maintenance and the ability to provide a solution for the broadest range of laboratory processes.



PERFORMANCE AND AFFORDABILITY - WITHOUT COMPROMISE

- Chemical resistant epoxy powder-coat
- Built with 16-gauge precision laser cut tubing
- Standard widths of 48", 60", 72", 84" and 96"
- 29"-37" height adjustable
- 3-tier single piece adjustable shelving made from 18 gauge steel
- Plumbing and electrical chase for easy installation and maintenance.

OPTIONAL EQUIPMENT

- 3 plumbing outlets per side
- Quick disconnect plumbing fittings
- Ceiling service panels with quick disconnect fittings
- Choice of leveling glides or casters
- Electrical receptacles
- Integrated lighting options
- Data receptacles
- Shared frame design for back-to-back bench configurations
- Rear and end modesty panel options
- Steel, epoxy, phenolic and laminate shelving options
- Front and rear seismic shelf lip options

